

### Introduction

The qualitative analysis, or identification, of the common anions is simpler than the analysis of the cations and usually depends on spot tests of the anions rather than separations followed by Confirmatory tests. Nevertheless, for the purpose of systematic qualitative analysis, anions are classified into 7 groups on the basis of their behaviour against  $\text{AgNO}_3$  and  $\text{BaCl}_2$  and the solubility in water and 2M  $\text{HNO}_3$  of precipitate products of the reactions.

**Group I:** anions form white precipitate of silver salts which are insoluble in water and 2M  $\text{HNO}_3$  solution. Barium salts are soluble in water. Group I consists of  $\text{Cl}^-$  (chloride),  $\text{Br}^-$  (bromide),  $\text{I}^-$  (iodide),  $[\text{Fe}(\text{CN})_6]^{4-}$  (ferrocyanide),  $[\text{Fe}(\text{CN})_6]^{3-}$  (ferricyanide) anions.

**Group II:** anions form white precipitates of Ag salts which are slightly soluble in water and soluble in 2M  $\text{HNO}_3$  solution. Barium salts are soluble in water. Group II consists of  $\text{NO}_2^-$  (nitrite),  $\text{CH}_3\text{COO}^-$  (acetate),  $\text{S}^{2-}$  (sulfide).

**Group III:** anions form white precipitates with both  $\text{Ag}^+$  and  $\text{Ba}^{2+}$  ions slightly soluble in water and soluble in 2M  $\text{HNO}_3$ . Group III consists of  $\text{SO}_3^{2-}$  (sulfite),  $\text{CO}_3^{2-}$  (carbonate),  $\text{C}_2\text{O}_4^{2-}$  (oxalate),  $\text{BO}_3^{3-}$  (borate),  $\text{C}_4\text{H}_4\text{O}_6^{2-}$  (tartrate) anions.

**Group IV:** anions form color precipitates with  $\text{Ag}^+$  and  $\text{Ba}^{2+}$  ions which are slightly soluble in water and soluble in 2M  $\text{HNO}_3$  solution. Group IV contain  $\text{PO}_4^{3-}$  (phosphate),  $\text{S}_2\text{O}_3^{2-}$  (thiosulfate),  $\text{CrO}_4^{2-}$  (chromate) anions.

**Group V:** anions do not form any precipitates with  $\text{Ag}^+$  and  $\text{Ba}^{2+}$  ions. This group consists of  $\text{NO}_3^-$  (nitric) and  $\text{MnO}_4^-$  (permanganate) anions.

**Group VI:**  $\text{SO}_4^{2-}$  ions form  $\text{Ag}_2\text{SO}_4$  precipitate soluble in water and white precipitate of  $\text{BaSO}_4$  insoluble in water and 2M  $\text{HNO}_3$ .

### Exercise 1:



#### Introduction

Group I anions form insoluble silver salts. Upon the addition of  $\text{AgNO}_3$   $\text{Cl}^-$ ,  $\text{Br}^-$ ,  $\text{I}^-$ ,  $[\text{Fe}(\text{CN})_6]^{4-}$ ,  $[\text{Fe}(\text{CN})_6]^{3-}$  ions precipitate as  $\text{AgCl}$ ,  $\text{AgBr}$ ,  $\text{AgI}$ ,  $\text{Ag}_4[\text{Fe}(\text{CN})_6]$  and  $\text{Ag}_3[\text{Fe}(\text{CN})_6]$  insoluble in 2 M  $\text{HNO}_3$

#### Separation and identification of I Group anions $\text{Cl}^-$ , $\text{Br}^-$ , $\text{I}^-$ , $[\text{Fe}(\text{CN})_6]^{4-}$ , $[\text{Fe}(\text{CN})_6]^{3-}$


The initial sample may be colorless or yellow. If the sample is yellow it may contain at least one of  $[\text{Fe}(\text{CN})_6]^{4-}$  and  $[\text{Fe}(\text{CN})_6]^{3-}$  ions. Therefore, they must be identified as first from the initial sample with the following reactions:

$[\text{Fe}(\text{CN})_6]^{4-}$




light blue color which changes into dark blue




$[\text{Fe}(\text{CN})_6]^{3-}$ 

reactions 3-4



reaction 5

After identification of these anions they have to be removed before further identification of the rest of anions of this group. For this purpose the excess of  $\text{ZnSO}_4$  solution should be used, because  $\text{ZnSO}_4$  reacts with  $[\text{Fe}(\text{CN})_6]^{4-}$  and  $[\text{Fe}(\text{CN})_6]^{3-}$  anions and precipitates as zinc salts of  $\text{Zn}_2[\text{Fe}(\text{CN})_6]$  and  $\text{Zn}_3[\text{Fe}(\text{CN})_6]$ :




reaction 6




reaction 7

When precipitates occurred, filtrate the mixture and discard the precipitate A (scheme 1). Filtrate A contains only  $\text{Cl}^-$ ,  $\text{Br}^-$  and  $\text{I}^-$  and should be colorless. If the Filtrate A is still yellow, additional portion of  $\text{ZnSO}_4$  must be added and resulted precipitate filtrated again. Then the rest of anions can be identified from obtained Filtrate A.

$\text{Br}^-$  and  $\text{I}^-$  ions can be also identified from the initial sample. The identification of these two ions can be followed in one test tube with the use of  $\text{Cl}_2\text{aq}$  and  $\text{CHCl}_3$ . For this purpose take 2 ml of the initial sample, add 2 ml of  $\text{CHCl}_3$  and some portion of  $\text{Cl}_2\text{aq}$  solution. Next, close test tube with your thumb and shake it vigorously. The presence of  $\text{I}^-$  ions is confirmed, when the organic layer of  $\text{CHCl}_3$  is tinted violet.




reaction 8



<http://www.slideshare.net/guest4f7e5c/redox-3462885>

The chlorine water should be added several time with small portions in order to check the presence of  $\text{Br}^-$ . When the violet color will disappear from the chloroform layer and yellow-orange color will appear we can state, that  $\text{Br}^-$  ions are present in our initial sample.




reaction 9

The lack of the color allow to state that  $\text{Br}^-$  ions are not present in the analyzed sample.

**NOTE:**  $\text{I}^-$  and  $[\text{Fe}(\text{CN})_6]^{3-}$  ions cannot be present together in one solution, because  $\text{I}^-$  ions undergo the oxidation in the presence of  $[\text{Fe}(\text{CN})_6]^{3-}$  ions.




Chloride ions  $\text{Cl}^-$  can be identified from initial sample. As first silver nitrate  $\text{AgNO}_3$  and 2 M  $\text{HNO}_3$  must be added. All anions from I group will precipitate as silver salts.



Silver chloride is insoluble in diluted inorganic acids and dissolves in diluted  $\text{NH}_3\text{aq}$  solution (1 portion of  $\text{NH}_3\text{aq}$  solution and 3 portions of distilled  $\text{H}_2\text{O}$ ):



After adding  $\text{HNO}_3$  solution to resulting Filtrate C (scheme 1) white precipitate of  $\text{AgCl}$  is obtained again:






### The rest of analytical reactions of $\text{Cl}^-$ , $\text{Br}^-$ , $\text{I}^-$ , $[\text{Fe}(\text{CN})_6]^{4-}$ , $[\text{Fe}(\text{CN})_6]^{3-}$ ions

#### $\text{H}_2\text{SO}_4$ diluted

$\text{Cl}^-$  no reaction

$\text{Br}^-$  and  $\text{I}^-$  react in higher temperature




#### $\text{KMnO}_4$ in acidic solution

$\text{Cl}^-$  no reaction in cold solution, but  $\text{Cl}^-$  undergo reduction after heating



$\text{Br}^-$ ,  $\text{I}^-$  and  $[\text{Fe}(\text{CN})_6]^{4-}$  undergo reduction at lower temperature




(and similarly for  $\text{I}^-$ )




#### **Cu(II) salts**

$\text{Cu(II)}$  salts undergo reduction in the presence of  $\text{I}^-$  ions and form  $\text{CuI}$  slightly soluble in water





**Hydrogen peroxide H<sub>2</sub>O<sub>2</sub>**


Iodide ions undergo oxidation in the presence of H<sub>2</sub>O<sub>2</sub> in the presence of weak acetic acid

**Fe(III) salts**





FeCl<sub>3</sub> evolve gaseous I<sub>2</sub> from I<sup>-</sup> ions solution



**Prussian blue** - a dark blue pigment

**Potassium dichromate K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub>**

causes the oxidation of Br<sup>-</sup> and I<sup>-</sup> in the presence of sulfuric acid H<sub>2</sub>SO<sub>4</sub>

**Other reducing agents in reactions with [Fe(CN)<sub>6</sub>]<sup>3-</sup> ions:**

Scheme 1

